

Short tutorial on ANACONDA, v2

- Installation
- Step-to-step through routines
- future merging with HDF5/MAX IV

Installation (as user)



- Download (clone) from GitHub:

https://github.com/gasfas/ANACONDA_2

- D&D in MATLAB cmd window

- Open 'GettingStarted.mlx',
or 'GettingStarted_plain.m'


- Run it, enjoy the errors...



Installation (as dev)

- **Install github desktop**
- **Clone from repository**
- **Request developer permissions**

Data analysis structure



Macro: function that executes subroutines.
One macro for each step: Correct, convert, plot, etc.

Importing (meta) data

Load data:

```
data = IO.import_raw('full data filename')
```

And metadata:

```
mdata = IO.import_metadata('full metadata filename')
```

If you have no data at hand, use an example:

```
data = IO.COBOLD.import_example();
```

If metadata is not yet defined: default metadata is defined in the toolbox for each setup. e.g.:

```
mdata = metadata.defaults.exp.CIEL.md_all_defaults
```

Quick visualization of imported data

Make histogram of X,Y:



plot.quickhist(data.h.det1.raw(:,1:2))

Make histogram of TOF:

plot.quickhist(data.h.det1.raw(:,3))


Inspect any function:

edit *macro.calibrate.TOF_2_m2q*



Correct data (t_0 , aberrations, etc.)

To handle alignment or artefacts issues



Correct data (t_0 , aberrations, etc.)

To handle alignment or artefacts issues

Calibration needed to determine the correction values (as defined in metadata)

For example to centre detector:


data_correct = macro.correct.dXdY(data, mdata, 'det1');

Or to deal with a spectrometer aberration:

data_correct = macro.correct.lensadd(data, mdata, 'det1')

On/off 'flags'

Corrections can be turned on or off in metadata:



```
/home/lna/Documents/Lund/ANACONDA_2_git/ANACONDA_2/package/+metadata/+defaults/+exp/+Laksman_TOF_e_XY/corr...
EDITOR PUBLISH VIEW
1 function [ exp_md ] = corr ( exp_md )
2 % This convenience function lists the default correction metadata, and can be
3 % read by other experiment-specific metadata files.
4 %% Correction parameters:
5 % Detector 1 %% Detector 1 %% Detector 1 %% Detector 1 %% Detector 1 %% Detector 1 %% Detector 1
6 % Which corrections should be performed on the data:
7 exp_md.corr.det1.ifdo.dxdy = true;% Does this data need detector image translation correction?
8 exp_md.corr.det1.ifdo.dtheta = true;% Does this data need detector image rotation correction?
9 exp_md.corr.det1.ifdo.dtof = true;% Does this data need detector absolute TOF correction?
10 exp_md.corr.det1.ifdo.R_circle = true;% is there non-roundness in the cylindrically symmetric data?
11
12 % The detector image translation parameters:
13 exp_md.corr.det1.dx = 0.0; %[mm] distance the center of detection is displaced left of the origin
14 exp_md.corr.det1.dy = 0.0; %[mm] distance the center of detection is displaced up the origin
15 % The detector image rotation parameters:
16 exp_md.corr.det1.dtheta = 0; %[deg] rotation of hits around the raw image centre (anticlockwise)
17 % The TOF deadtime correction parameter:
18 exp_md.corr.det1.dtof = 0;% [ns] The difference between signal propagation times of trigger and hit
19
20 exp_md.corr.det1.R_circle = load(fullfile(fileparts(mfilename('fullpath'))), 'R_circle_param.mat');
21
22 % Detector 2 %% Detector 2 %% Detector 2 %% Detector 2 %% Detector 2 %% Detector 2 %% Detector 2
23
24 % Which corrections should be performed on the data:
25 exp_md.corr.det2.ifdo.dxdy = true;% Does this data need detector image translation correction?
26 exp_md.corr.det2.ifdo.dtheta = true;% Does this data need detector image rotation correction?
```

On/off 'ifdo' flags

These flags exist for the macros:

- Correct
- Convert
- Plot
- Filter

Convert (X, Y, TOF, etc.) to (m, p_x , etc.)

```
/home/ina/Documents/Lund/ANACONDA_2_git/ANACONDA_2/package/+metadata/+defaults/+exp/+Laksman_TOF_e_XY/conv.m
EDITOR PUBLISH VIEW
1 function [ exp_md ] = conv ( exp_md )
2 % This convenience function lists the default conversion metadata, and can be
3 % read by other experiment-specific metadata files.
4
5 %% Conversion factors:
6 % Which conversions should be performed on the data:
7
8 exp_md.conv.det2.ifdo.m2q = true;% Does the user want to convert to mass-over-charge?
9 exp_md.conv.det2.ifdo.m2q_label = true;% Does the user want to convert to mass-2-charge labels
10 exp_md.conv.det2.ifdo.m2q_group_label = true;% Does the user want to convert to mass-2-charge groups
11 exp_md.conv.det2.ifdo.cluster_size = false;% Does the user want to convert to cluster size?
12 exp_md.conv.det2.ifdo.momentum = true;% Does the user want to convert to momentum?
13 exp_md.conv.det2.ifdo.KER = true;% Does the user want to convert to Kinetic energy?
14 exp_md.conv.det2.ifdo.R_theta = true;% Does the user want to convert to R-theta coordinates?
15 exp_md.conv.det2.ifdo.angle_p_corr_C2 = true;% Does the user want to calculate mutual momentum angles
16 exp_md.conv.det2.ifdo.angle_p_corr_C3 = true;% Does the user want to calculate mutual momentum angles
17 exp_md.conv.det2.ifdo.q_label = false;% Does the user want to convert to have q-labels define
18 exp_md.conv.det2.ifdo.CSD = false;% Does the user want to convert to CSD?
19 % Detector 1 %% Detector 1 %% Detector 1 %% Detector 1 %% Detector 1 %% Detector 1 %% Detector 1
20
21 % Time Of Flight to Mass to Charge calibration
22 exp_md.conv.det2.m2q.factor = 962.3;%Conversion factor from mass 2 charge to TOF
23 exp_md.conv.det2.m2q.t0 = 20.5;%[ns] dead time correction
24
25 % Mass to charge to Mass to Charge labels
Ln 1 Col 1
```

Filter on hits:

Filters are formed from *conditions*. e.g.:

```
cond.oil.data_pointer = 'h.det1.m2q_l';
```

Filter on hits:

Filters are formed from *conditions*. e.g.:

```
cond.oil.data_pointer      = 'h.det1.m2q_l';  
cond.oil.type              = 'discrete';
```

Filter on hits:

Filters are formed from *conditions*. e.g.:

```
cond.oil.data_pointer    = 'h.det1.m2q_l';  
cond.oil.type            = 'discrete';  
cond.oil.value           = [72; 73];
```

Filter on hits:

Filters are formed from *conditions*. e.g.:

cond.oil.data_pointer	= 'h.det1.m2q_l';
cond.oil.type	= 'discrete';
cond.oil.value	= [72; 73];
cond.oil.translate_condition	= 'AND';

Filter on hits:

Filters are formed from *conditions*. e.g.:

cond.oil.data_pointer	= 'h.det1.m2q_l';
cond.oil.type	= 'discrete';
cond.oil.value	= [72; 73];
cond.oil.translate_condition	= 'AND';
cond.oil.invert_filter	= true;

Filter on hits:

Data pointer could be any MATLAB expression:

```
cond.m2q.data_pointer    = 'h.det1.m2q_l.*h.det1.TOF+1';  
cond.m2q.type            = 'discrete';  
cond.m2q.value           = [72; 73];  
cond.m2q.translate_condition = 'AND';  
cond.m2q.invert_filter   = true;
```

Filter on hits:

The 1st and/or 2nd and/or 3rd hit can be selected:

```
cond.m2q.data_pointer    = 'h.det1.m2q_l.*h.det1.TOF+1';  
cond.m2q.type            = 'discrete';  
cond.m2q.value           = [72; 73];  
cond.m2q.translate_condition = 'AND';  
cond.m2q.invert_filter   = true;  
cond.m2q.hitselect       = [1, 3];
```

Filter on events:

Total Kinetic energy:

```
KER_sum.type          = 'continuous';  
KER_sum.data_pointer = 'e.det1.KER_sum';  
KER_sum.value         = [0; 80]; % [eV]
```

Merging of conditions:

When merging filters, operators can be defined:


```
cond.combined.m2q      = m2q  
cond.combined.oil      = oil;  
cond.combined.KER_sum  = KER_sum;  
cond.combined.operators = {'AND', 'OR'};
```

Making Plots

Metadata separated in categories:

figure: Figure handle properties
axes: Axes handle properties
hist: Histogram metadata
GraphObj: Graphical Object handle
Cond: (optional) plot condition


What are handles?

A screenshot of the MATLAB Command Window. The title bar reads "Command Window". The command prompt shows the command `>> mdata.plot.det1.m2q` being executed. The output is `ans =` followed by a `struct` with fields: `figure: [1x1 struct]`, `axes: [1x2 struct]`, `hist: [1x1 struct]`, and `GraphObj: [1x1 struct]`. The cursor is at the prompt `fx >> |`.

```
Command Window
>> mdata.plot.det1.m2q
ans =
    struct with fields:
        figure: [1x1 struct]
         axes: [1x2 struct]
         hist: [1x1 struct]
    GraphObj: [1x1 struct]
fx >> |
```

Calibration

All correction and conversion parameters need calibration.




Calibration

Metadata calibration:

`macro.calibrate.momentum(data, 1, mdata)`

HDF5 dummy file



HDFGridView: events				HDFGridView: hits			
Edit				Edit			
	0	1	2		0	1	2
0	1	16374549200	0	0	14.2357106652	-15.4965712688	0.0
1	2	33498408400	1	1	-16.4575149984	27.3239613508	0.0
2	3	36610398800	2	2	25.0500673308	7.32857051426	0.0
3	4	59108469200	3	3	-9.79397966569	15.1760362528	0.0
4	5	66322050000	4	4	32.8930096634	9.36751911842	0.0
5	6	75185074000	5	5	-3.31633399967	9.43021992748	0.0
6	7	111851582800	6	6	0.075369333258	0.628233227764	0.0
7	8	124355662800	7	7	2.53150399975	-0.448619580562	0.0
8	9	127945448400	8	8	-24.4167613309	16.6469022824	0.0
9	10	171710728400	9	9	-6.00801666607	-3.63881314366	0.0
10	11	177175048400	10	10	-37.3707653296	5.99655316627	0.0
11	12	187283429200	11	11	-1.24335199988	20.8794635264	0.0
12	13	197100491600	12	12	-8.7081016658	32.1184159752	0.0
13	14	205114933200	13	13	-34.4201596632	-15.870497322	0.0
14	15	212016181200	14	14	22.0428933311	20.6355353763	0.0
15	16	226913973200	15	15	27.5718726639	9.11326621472	0.0
16	17	227703774800	16	16	24.8353279975	18.8851342777	0.0
17	18	228904549200	17	17	-22.7428959977	12.1271300985	0.0
18	19	235725400400	18	18	8.06141733253	-4.75173763778	0.0
19	20	246668731600	19	19	-29.2130559971	13.5742699217	0.0
20	21	291957614800	20	20	26.8788703306	20.2742174968	0.0
21	22	312806562000	21	21	-34.6499339965	-3.2473743269	0.0
22	23	351951115600	22	22	6.51179566602	26.8627670829	0.0
23	24	358964443600	23	23	28.4108796638	-7.21351097026	0.0
24	25	371575765200	24	24	-10.209241999	24.7262613061	0.0
25	26	374717800400	25	25	20.6262159979	22.4102821781	0.0
26	27	397279154000	26	26	23.3811866643	-25.7729360194	0.0
27	28	410704558800	27	27	22.5126859977	-25.8592602666	147.938
28	29	413822053200	28	28	-23.1197886644	11.7732248183	0.0
29	30	434146174800	29	29	12.3361506654	-3.3943642237	0.0
				30	-31.8734816635	2.03217744101	0.0